



NEW 3D PRINTED JOEYS

Joeys



PRESCRIPTION TOOL
Summary of Features

Joey Type	Heel Type	Correction	Arch Fill Technique	Included Cast & Shell Modifications
Court Hook	Hook	4°	Mid	
Support	Concave wedge	4°	Mid	
Support Plus	Concave wedge	4°	Mid	Approx. 3mm higher arch than 'Joeys Support'
Lateral Support	Concave wedge	5° Valgus	Mid/Mod Root	Plantar 5th ray grind Cuboid notch
Original	Wedge (Now standard with heel stabiliser)	13°	Inverted/ Mid	1st ray accommodation Cuboid notch
Original Plus	Wedge (Now standard with heel stabiliser)	13°	Inverted/ Mid	1st ray accommodation Cuboid notch Medial flare
Inverted	Wedge (Now standard with heel stabiliser)	20°	Inverted	1st ray accommodation Cuboid notch Medial flare
Junior	Traditional with high medial lateral heel cup	0°	Mid	Medial wrap Long lateral heel cup
Plantar Heel Pain	Traditional	2°	Mid/Mod Root	Sweet spot (medial tubercle) Plantar fascial accommodation 1st ray accommodation Plantar 5th ray grind

NEED MORE INFORMATION?

For more information about our brand new range of 3D printed Joeys, please call; **1300 667 744** or email us direct at **info@footwork.com.au**.

www.footwork.com.au